

KERN 362-990-300

KERN

Wooden box for weight sets (Cylindrical/compact, 1 g - 10 kg)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	M3
Case suitable for nominal weight values	1 g - 10 kg
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	510×150×70 mm
Material	beech
Case - Smallest possible weight	1 g
Case - Highest possible weight	10 kg

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	510×150×70 mm
Net weight	2,5 kg
Shipping method	Parcel service
Net weight approx.	2,6 kg
Gross weight approx.	2,6 kg
Shipping weight	2,5 kg

Pictograms

STANDARD

