

Product brochure

07.08.2026

KERN 366-86

Individual weight, OIML Class M3, 5 kg, Block, cast iron lacquered

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
317-280	Gloves 317-280 (cotton, 1 pair)
346-060-600	Aluminium weight case 346-060-600 (for nominal values 5 kg Block for classes F1 - M3 for design Block)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-643	DAkkS calibration certificate 962-643
962-643R	DAkkS Recalibration 962-643R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days