

Product brochure

07.08.2026

KERN 366-89

Individual weight, OIML Class M3, 50 kg, Block, cast iron lacquered

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

PICTOGRAMS

Standard

Option





ACCESSORIES

Model	Description
317-290	Gloves 317-290 (leather, cotton, 1 pair)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-646	DAkkS calibration certificate 962-646
962-646R	DAkkS Recalibration 962-646R

ENGLISH

KERN Pictogram Overview



Pallet shipment: The duration of internal product provision are two days



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days