

KERN 366-91

KERN

Individual weight, OIML Class M3, 100 g, Cylinders, cast iron lacquered



Dimensions packaging (W×D×H)	31×31×22 mm
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg

Pictograms

STANDARD



Category	
Brand	KERN
Product category	Individual weight
Product group	OIML M3
Product family	366

Measuring System	
Conformity assessment possible	no
Test weight - volume (cm ³)	14,0845
Test weight - volume uncertainty (cm ³)	1,1902
Test weight - Density (kg/m ³)	7100
Test weight - Density uncertainty (kg/m ³)	600
Test weight - Density [Min] (OIML)	1500
Test weight - nominal value	100 g
OIML-compliant	✓
Weight OIML class	M3
Permissible tolerance according to OIML (+/-)	50 mg

Construction	
Shape weight	Cylinders
Dimensions (Ø×H)	31×22 mm
Material	cast iron lacquered

Services (optional)	
Article number for DAkkS calibration	962-637

Packing & Shipping	
Delivery time	1 d