

KERN 366-95

Individual weight, OIML Class M3, 2 kg, Knob, cast iron lacquered



Dimensions packaging (W×D×H)	70×70×113 mm
Shipping method	Parcel service
Net weight approx.	2,0 kg
Gross weight approx.	2,0 kg
Shipping weight	2 kg

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M3
Product family	366

Measuring System

Conformity assessment possible	no
Test weight - volume (cm ³)	281,6901
Test weight - volume uncertainty (cm ³)	23,8048
Test weight - Density (kg/m ³)	7100
Test weight - Density uncertainty (kg/m ³)	600
Test weight - Density [Min] (OIML)	1500
Test weight - nominal value	2 kg
OIML-compliant	✓
Weight OIML class	M3
Permissible tolerance according to OIML (+/-)	1 g

Construction

Shape weight	Knob
Dimensions (Ø×H)	70×113 mm
Material	cast iron lacquered

Services (optional)

Article number for DAkkS calibration	962-642
--------------------------------------	---------

Packing & Shipping

Delivery time	1 d
---------------	-----