

KERN ATB-US03



Ultrasound contact gel, refill pack, approx. 70 ml



Category

| | |
|------------------|-----------------------|
| Brand | Sauter |
| Product category | Measurementtechnology |

Packing & Shipping

| | |
|------------------------------|----------------|
| Delivery time | 1 d |
| Dimensions packaging (W×D×H) | 150×100×2 mm |
| Net weight | 0,016 kg |
| Shipping method | Parcel service |
| Net weight approx. | 0,05 kg |
| Gross weight approx. | 0,10 kg |
| Shipping weight | 0,1 kg |

Pictograms

STANDARD

