


# SAUTER CE Q24



Traction device for "S" measuring cell M24×2, safe load of 8150 kg



**Date of creation: 20.06.2026**

## Category

Brand	Sauter
Product category	Measurementtechnology
Product group	Traction device

## Construction

Material	steel, galvanized
----------	-------------------

## Packing & Shipping

Delivery time	1 d
Shipping method	Parcel service
Net weight approx.	8 kg
Gross weight approx.	8 kg
Shipping weight	7,805 kg

## Product information

GTIN/EAN number	4045761221697
Trade number	84239010
REACH classification	no prohibited substances

## Pictograms

### STANDARD

