

SAUTER CE Q30906

Elastomer bearing with hold-down device

KERN



20.06.2026

Category

| | |
|------------------|-----------------------|
| Brand | Sauter |
| Product category | Measurementtechnology |
| Product group | Elastomer bearing |

Construction

| | |
|----------|------------------|
| Material | steel, stainless |
|----------|------------------|

Packing & Shipping

| | |
|----------------------|----------------|
| Delivery time | 1 d |
| Shipping method | Parcel service |
| Net weight approx. | 0,50 kg |
| Gross weight approx. | 0,55 kg |
| Shipping weight | 0,502 kg |

Product information

| | |
|----------------------|--------------------------|
| GTIN/EAN number | 4045761221758 |
| Trade number | 84239010 |
| REACH classification | no prohibited substances |

Pictograms

STANDARD

