


SAUTER CE Q34903

Base plate. Galvanised steel

KERN



Date of creation: 20.06.2026

Category

Brand	Sauter
Product category	Measurementtechnology
Product group	Plate

Construction

Material	steel, galvanized
----------	-------------------

Packing & Shipping

Delivery time	1 d
Shipping method	Parcel service
Net weight approx.	3,0 kg

Product information

GTIN/EAN number	4045761221789
Trade number	84239010
REACH classification	no prohibited substances

Pictograms

STANDARD

