

# KERN CKE 36K0.1-Inventory

Carry out stocktaking quickly and accurately using a combination of a counting balance and software, also ideal for large warehouses with vast numbers of parts



Software EasyTouch

SET-52  
Inventory

Smallest piece weight when piece counting - normal conditions 1 g

## Interface

EasyTouch (ET) - Compatibility Vollständig

## Power Supply

Battery operating time 20 h  
 Rechargeable battery operating time - backlight on 24 h  
 Rechargeable battery operating time - backlight off 48 h  
 Rechargeable battery charging time 8 h

## Environmental conditions

Ambient temperature [Min] -10 °C  
 Ambient temperature [Max] 40 °C  
 Humidity of environment [Max] 80 %

## Approval

CE mark ✓

## Services (optional)

Article number for DAkKS calibration 963-128  
 Article number for certificate of conformity (verification) 969-517

## Packing & Shipping

Delivery time 1 d  
 Dimensions packaging (W×D×H) 470×470×190 mm  
 Shipping method Parcel service  
 Net weight approx. 7 kg  
 Gross weight approx. 8 kg  
 Shipping weight 8,4 kg

## Category

Brand KERN  
 Product category Balance software bundle  
 Product group Stocktake Scale  
 Product family CKE-Inventory

## Measuring System

Weighing system Strain gauge  
 Weighing capacity [Max] 36 kg  
 Readability [d] 100 mg  
 Reproducibility 200 mg  
 Linearity ± 500 mg  
 Resolution 360.000  
 Recommended adjusting weight 30 kg (E2)  
 Possible calibration points 10 kg; 20 kg; 30 kg  
 Stabilization time 3 s  
 Warm-up time 120 min  
 Eccentric loading at 1/3 [Max] 1 g  
 Maximum creep (15 minutes) 1 g  
 Maximum creep (30 minutes) 2 g

## Construction

Dimension housing (W×D×H) 350×390×120 mm  
 Dimensions weighing surface (W×D) 340×240 mm  
 Dimensions weighing surface 340×240 mm  
 Dimensions weighing platform (W×D×H) 340×240×21 mm  
 Level indicator ✓  
 Levelling feet adjustable ✓

## Functions

Counting resolution (laboratory conditions) 360.000  
 Counting reference weight can be entered ✓  
 Smallest piece weight when piece counting - laboratory conditions 100 mg

# KERN CKE 36K0.1-Inventory



Carry out stocktaking quickly and accurately using a combination of a counting balance and software, also ideal for large warehouses with vast numbers of parts

## Pictograms

### STANDARD



### OPTION

