

KERN DS-A01

KERN

Set for underfloor weighing, consists of platform, bow, hook, only for models with weighing plate
308×318×75 mm



Category

| | |
|------------------|---------------------|
| Brand | KERN |
| Product category | Balance |
| Product group | Underfloor weighing |

Construction

| | |
|----------|-----------------|
| Material | Stainless steel |
|----------|-----------------|

Packing & Shipping

| | |
|------------------------------|----------------|
| Delivery time | 1 d |
| Dimensions packaging (W×D×H) | 340×330×60 mm |
| Net weight | 2 kg |
| Shipping method | Parcel service |
| Net weight approx. | 2,0 kg |
| Gross weight approx. | 2,0 kg |
| Shipping weight | 2 kg |

Pictograms

STANDARD

