


KERN FK-A02

KERN

Tensiometer kit for high-capacity tensile strength testing up to 1000 N, steel attachment and steel rolls, rolls cannot be adjusted



Category

| | |
|------------------|-----------------------|
| Brand | Sauter |
| Product category | Measurementtechnology |

Construction

| | |
|---------------------------------------|---------------|
| Dimensions device (W×D×H) | 140×140×40 mm |
| Dimensions completely mounted (W×D×H) | 152×80×28 mm |
| Material | metal |
| Material housing | metal |

Environmental conditions

| | |
|---------------------------|--------|
| Storage temperature [Min] | -10 °C |
| Storage temperature [Max] | 40 °C |

Packing & Shipping

| | |
|----------------------|----------------|
| Delivery time | 1 d |
| Net weight | 0,13 kg |
| Shipping method | Parcel service |
| Net weight approx. | 0,15 kg |
| Gross weight approx. | 0,15 kg |
| Shipping weight | 0,133 kg |

Pictograms

STANDARD

