

KERN HO-A08

KERN

Test stand for repeatable movements during testing



Category

| | |
|------------------|-----------------------|
| Brand | Sauter |
| Product category | Measurementtechnology |

Construction

| | |
|---------------------------|----------------|
| Dimension housing (W×D×H) | 200×200×332 mm |
| Material | metal |
| Material housing | metal |

Packing & Shipping

| | |
|------------------------------|----------------|
| Delivery time | 1 d |
| Dimensions packaging (W×D×H) | 320×320×440 mm |
| Net weight | 8,8 kg |
| Shipping method | Parcel service |
| Net weight approx. | 9 kg |
| Gross weight approx. | 9 kg |
| Shipping weight | 9 kg |

Pictograms

STANDARD

