

KERN MAP 80K1S

KERN

Hand grip dynamometer, e.g. for rehabilitation treatment after accidents



Category	
Brand	KERN
Product category	Medical scale
Product group	Hand grip dynamometer
Product family	MAP

Measuring System	
Weighing capacity [Max]	80 kg
Readability [d]	100 g
Recommended adjusting weight	40 kg (M1)
Possible calibration points	40 kg
Units	kg lb

Display	
Display type	LCD graphic
Display digit height	12 mm

Construction	
Dimensions device (W×D×H)	55×88×212 mm

Functions	
Auto-off interval(s) in battery mode/rechargeable battery mode	1 min

Power Supply	
Supplied power supply	Battery
Battery / accumulator type	Li-Metal
Battery	CR2450
Battery operating time	53 h
Battery capacity	630 mAh
Battery voltage	3 V

Environmental conditions	
Humidity of environment [Max]	80 %
Use temperature [Min]	5 °C
Use temperature [Max]	35 °C
Storage temperature [Min]	-20 °C
Storage temperature [Max]	60 °C

Services (optional)	
Article number for factory calibration	961-102K

Packing & Shipping	
Dimensions transport container (W×D×H)	350×265×85 mm
Delivery time	1 d
Dimensions packaging (W×D×H)	367×285×100 mm
Shipping method	Parcel service
Net weight approx.	0,35 kg
Gross weight approx.	1,8 kg
Shipping weight	2,1 kg

Pictograms

STANDARD



OPTION

