

KERN MBC-A05

KERN

Foot/head rest to fix to the lying surface for even more safety when weighing small children, scope of delivery: 1 head rest and 1 foot rest



Pictograms

STANDARD



Date of creation: 19.06.2026

Category

Brand	KERN
Product category	Balance
Product group	Foot/head rest

Construction

Dimension (W×D×H)	600×280×88 mm
Dimensions completely mounted (W×D×H)	600×280×128 mm
Dimensions weighing surface (W×D)	535×280 mm

Environmental conditions

Storage temperature [Min]	-20 °C
Storage temperature [Max]	60 °C

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	650×300×155 mm
Shipping method	Parcel service
Net weight approx.	1,4 kg
Gross weight approx.	2,0 kg
Shipping weight	6 kg
Shipping volume weight	6 kg

Product information

GTIN/EAN number	4045761158610
Trade number	84239010
REACH classification	no prohibited substances