

# KERN MCD 300K-1

**KERN**

Mobile chair scale in ergonomically-optimised design for safe, convenient weighing



## Category

Brand	KERN
Product category	Medical scale
Product group	Chair scale
Product family	MCD

## Measuring System

Construction type of the scale	Single-range balance
Weighing system	Strain gauge
Weighing capacity [Max]	300 kg
Readability [d]	0,1 kg
Reproducibility	0,1 kg
Linearity	± 0,1 kg
Resolution	3.000
Tare range	299,9 kg
Adjustment options	Adjusting with external weight
Recommended adjusting weight	300 kg (M1)
Stabilization time	3 s
Warm-up time	10 min
Eccentric loading at 1/3 [Max]	0,1 kg
Maximum creep (15 minutes)	400 g
Maximum creep (30 minutes)	800 g
Default unit	kg
Units	kg

## Display

Display type	LCD
Display digit height	25 mm
Languages of the user interface	Symbol language

## Construction

Dimensions completely mounted (W×D×H)	625×990×985 mm
Dimension display device (W×D×H)	200×130×60 mm
Dimensions weighing surface (W×D)	390×360 mm
Dimensions weighing surface	390×360 mm
Dimensions weighing platform (W×D×H)	500×360×370 mm
Material housing	steel, lacquered
Material weighing plate	plastic
Material platform	Plastics
Material display housing	Plastics
Level indicator	✓
Transport wheels	✓

## Functions

Function for BMI calculation	✓
Hold function	✓
Auto-off interval(s) in battery mode/rechargeable battery mode	off 3, 5, 15, 30 min
Tare function	manual (single & multi)
Number of keys for operation	6

## Power Supply

Supplied power supply	Battery
Battery / accumulator type	Alkali(-Manganese)
Battery	6×1.5 V AA
Battery operating time	40 h
Battery capacity	55 mAh
Battery voltage	18 V

## Environmental conditions

Ambient temperature [Min]	10 °C
Ambient temperature [Max]	40 °C
Humidity of environment [Max]	80 %

## Approval

CE mark	✓
---------	---

## Services (optional)

Article number for adjustment at the location of installation	961-249
Article number for DAkkS calibration	963-129
Article number for certificate of conformity (verification)	969-517

# KERN MCD 300K-1



Mobile chair scale in ergonomically-optimised design for safe, convenient weighing

## Packing & Shipping

Delivery time	2 d
Dimensions packaging (W×D×H)	1070×690×1080 mm
Net weight	22 kg
Shipping method	Freight
Net weight approx.	22 kg
Gross weight approx.	36 kg
Shipping weight	35 kg

## Pictograms

### STANDARD



### OPTION

