



The robust classic for quick determination of body weight

Features

- Timeless design
- Robust construction
- **1** Large, clear circular scale which is easy to read
- **2** Large tread with anti-slip surface made of plastic
- Secure and non-slip positioning with rubber feet
- Easy to operate: stand – read the weight
- Easy adjustment of the display using the setting wheel
- Delivered in single design packaging
- **Note:** The model with the suffix -S05 is delivered in a set of 5 units. i.e. the price given in the table refers to 5 items

Technical data

- Digit height 8 mm
- Overall dimensions W×D×H 290×425×90 mm
- Net weight approx. 3,2 kg

Accessories

- Alcohol-free cleaning cloths for disinfectant wiping, quick acting, based on modern quaternary ammonium compounds and effective against papova viruses. Particularly gentle on materials. Fulfill the legal requirements for occupational safety in accordance with TRGS 525/540. Packaging contents 100 pcs., KERN MYC-01

STANDARD



Model	Weighing capacity	Readability	Quantity delivered (balance)
	[Max] kg	[d] kg	
KERN			
MGC 100K-1	150	1	1
MGC 100K-1S05	150	1	5