



Discover more details and matching accessories online!

The robust classic for quick determination of body weight

Features

- Timeless design
- Robust construction
- **1** Large, clear circular scale which is easy to read
- **2** Large tread with anti-slip surface made of plastic
- Secure and non-slip positioning with rubber feet
- Easy to operate: stand - read the weight
- Easy adjustment of the display using the setting wheel
- Delivered in single design packaging

Technical data

- Digit height 8 mm
- Overall dimensions W×D×H 290×425×90 mm
- Net weight approx. 3,2 kg

Note: The model with the suffix -S05 is delivered in a set of 5 units. i.e. the price given in the table refers to 5 items

STANDARD



Model	Weighing capacity [Max] kg	Readability [d] kg	Quantity delivered (balance)
KERN			
MGC 100K-1	150	1	1
MGC 100K-1S05	150	1	5