



Discover more details and matching accessories online!

Design body analysis scale to determine body weight, body fat, body fluid, muscle mass, proportion of bone, BMI and metabolic rate/metabolic rate when exercising

Features

- Used to determine body fat, body water, muscle mass with an accuracy of 0.1 % as well as the body weight and bone mass in 0.1 kg increments
- Display showing the calorie requirements for maintaining the current body weight depending on the level of activity and sex
- Sturdy glass weighing plate with four large contact surfaces for maximum precision
- Easy and hygienic cleaning
- Particularly flat design facilitates access
- Secure and non-slip positioning with rubber feet
- Easy to operate: simply tap on it – stand on it – read the weight
- Memory for up to 10 people (gender, age, height, level of activity)
- Integrated AUTO-OFF function to preserve the batteries
- Weight displayed in kg or lb
- **1** Large high-contrasted display that is easy to read
- **2** Delivered in single design packaging

Technical data

- LCD display, digit height 35 mm
- Overall dimensions W×D×H 300×300×25 mm
- Batteries included, 2×CR2032, operating time up to 100 h
- Net weight approx. 1,8 kg

Note: The model with the suffix -S05 is delivered in a set of 5 units. i.e. the price given in the table refers to 5 items

STANDARD



Model	Weighing capacity	Readability	Quantity delivered (balance)
	[Max] kg	[d] kg	
KERN			
MFB 150K100	182	0,1	1
MFB 150K100S05	182	0,1	5