



Discover more details and matching accessories online!

Personal scale with particularly large, flat weighing plate to make it easier to step on and give maximum safety

Features

- This personal scale is the ideal measuring instrument for elderly care, paramedics or private sector
- Very large weighing plate
- Easy and hygienic cleaning
- **1** The scale can be easily transported using a handle and does not require much storage space
- **2** Very large platform with non-slip and wear-resistant plastic-surface
- Particularly flat design facilitates access
- Easy to operate: simply tap on it – stand on it – read the weight
- Integrated AUTO-OFF function to preserve the batteries
- Weight displayed in kg or lb
- Large display that is easy to read
- Delivered in single design packaging

Technical data

- Large LCD display, digit height 30 mm
- Overall dimensions WxDxH 434x370x40 mm
- Ready for use: Batteries included, 4x1,5 V AAA, operating time up to 20 h
- Net weight approx. 3,2 kg

Note: The model with the suffix -S05 is delivered in a set of 5 units. i.e. the price given in the table refers to 5 items

STANDARD



Model	Weighing capacity		Readability		Quantity delivered (balance)
	[Max]	kg	[d]	kg	
KERN					
MPI 200K-1	200		0,1		1
MPI 200K-1S05	200		0,1		5