



Discover more details and matching accessories online!

Personal scale with column for convenient, ergonomically-optimised operation

Features

- This personal scale is the ideal measuring instrument for elderly care, paramedics, fitness clubs or in the private sector and can weigh overweight people up to 250 kg
- KERN MPL: BMI function to determine underweight/normal weight/surplus weight
- KERN MPL: Tare function to zero the balance
- Weight displayed in kg or lb
- **1** Particularly flat design facilitates access
- Sturdy glass weighing plate for secure, comfortable standing position
- Easy and hygienic cleaning
- **2** Large display that is easy to read
- **2** The column, which is firmly fixed to the scale, makes it easy to operate the scale as well as to read the weighing results even under poor lighting conditions
- Easy to operate:
simply tap on it – stand on it – read the weight
- Integrated AUTO-OFF function to preserve the batteries
- **3** Delivered in design packaging

Technical data

- Large LCD display, digit height 28 mm
- Weighing platform dimensions W×D×H 350×350×45 mm
- Overall dimensions W×D×H 350×420×955 mm (incl. stand)
- Batteries included, 4×1.5 V AA, operating time up to 100 h
- Net weight approx. 3,8 kg

STANDARD



Model	Weighing capacity	Readability
	[Max] kg	[d] kg
KERN		
MPK 200K-1P	250	0,1
MPL 200K-1P	250	0,1