

KERN MSC 100

KERN

When every centimetre counts – mechanical height rod approved as a medical device for professional use in medical diagnostics



Pictograms

STANDARD



Category

Brand	KERN
Product category	Medical scale
Product group	Height rod
Product family	MSC

Measuring System

Measuring range measuring rod (cm)	7,5-100cm
Units	inch/ft cm
Measuring range length [Min]	7,5 cm
Measuring range length [Max]	100 cm
Tolerance measuring rod	5 mm
Division	1 mm

Construction

Dimension housing (W×D×H)	1020×290×50 mm
Dimensions device (W×D×H)	1040×290×55 mm

Environmental conditions

Use temperature [Min]	10 °C
Use temperature [Max]	40 °C
Storage temperature [Min]	-20 °C
Storage temperature [Max]	60 °C

Approval

CE mark	✓
---------	---

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	1085×70×40 mm
Net weight	0,694 kg
Shipping method	Parcel service
Net weight approx.	0,70 kg
Gross weight approx.	0,95 kg
Shipping weight	0,919 kg