



Discover more details and matching accessories online!

Convenient circumference measuring tape with practical pull-out mechanism

Features

- Continuous pull-out measuring tape for precise measurement e.g. around the body
- Ideal to determine the WHR (Waist-to-Hip-Ratio), which can give information about the distribution of the body's fat deposits using an indirect determination of abdominal fat
- With integrated feed mechanism
- The easiest handling:
 - Pull out the measuring tape
 - Hang it in the recess provided for this
 - Operate the feed mechanism
 - Read measurement value with millimetre precision

Technical data

- Measuring tape length 205 cm (80 inch)

Note: The model with the suffix -S05 is delivered in a set of 5 units. i.e. the price given in the table refers to 5 items

STANDARD



Model	Measuring range	Increments of	Overall dimensions	Net weight
	mm (inch)	mm (inch)	W×D×H mm	approx. kg
KERN MSW 200	10 – 205 (4 – 80) inch	1 (0,1)	55×90×20	0,05
MSW 200S05	10 – 205 (4 – 80) inch	1 (0,1)	55×90×20	0,25