

# KERN OBB-A 1 103

**KERN**

Abbe N.A. condenser 1,25 precentred (with aperture diaphragm)



## Category

Brand	Optics
Product category	Optical Instrument
Product group	Microscope condenser

## Condensor

Condensor type	Abbe
----------------	------

## Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	80×80×70 mm
Shipping method	Parcel service
Net weight approx.	0,30 kg
Gross weight approx.	0,80 kg
Shipping weight	0,8 kg

## Product information

GTIN/EAN number	4045761164420
Trade number	90119000
REACH classification	no prohibited substances

## Pictograms

STANDARD

