

KERN OBB-A 1242

KERN

Planachromatic objective lens 100 x / 1,25 (0,69 mm) (sprung) (oil)



Pictograms

STANDARD



Category

Brand	Optics
Product category	Optical Instrument
Product group	Microscope objective lens

Construction

Anti fungus	✓
-------------	---

Objective

Objective magnification	100 x
Objective quality	Plan achromatic
Numerical Aperture	1,25
Objective working distance	0,69 mm

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	65×40×40 mm
Net weight	0,3 kg
Shipping method	Parcel service
Net weight approx.	0,30 kg
Gross weight approx.	0,60 kg
Shipping weight	0,6 kg