

KERN OBB-A 1479

KERN

Achromatic objective lens, 60 x /0,85 (sprung) W.D. (0,1 mm)



Pictograms

STANDARD



Category

Brand	Optics
Product category	Optical Instrument
Product group	Microscope objective achromatic

Construction

Anti fungus	✓
-------------	---

Objective

Objective magnification	60 x
Objective quality	Achromatic
Numerical Aperture	0,85
Objective working distance	0,1 mm

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	65×40×40 mm
Net weight	0,3 kg
Shipping method	Parcel service
Net weight approx.	0,30 kg
Gross weight approx.	0,60 kg
Shipping weight	0,6 kg