

# KERN OBB-A 1503

Specimen holder (Ø 110 mm)

**KERN**



## Category

Brand	Optics
Product category	Optical Instrument

## Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	130×130×10 mm
Net weight	0,1 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,20 kg
Shipping weight	0,2 kg

## Pictograms

STANDARD

