

KERN OBB-A 1589

KERN

LED spare bulb



Date of creation: 22.06.2026

Category

Brand	Optics
Product category	Optical Instrument
Product group	Microscope illumination

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	70×25×25 mm
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,20 kg
Shipping weight	0,2 kg
Shipping volume weight	0,1 kg

Product information

GTIN/EAN number	404576 1249332
Trade number	85395200
REACH classification	no prohibited substances

Pictograms

STANDARD

