

KERN OBB-A 1652

KERN

5W LED Fluorescence unit (B / G) for inverted microscope (KERN OCM 167)



Pictograms

STANDARD



Date of creation: 23.06.2026

Category

Brand	Optics
Product category	Optical Instrument
Product group	Fluorescence unit

Construction

Material	metal, lacquered
----------	------------------

Power Supply

Supplied power supply	Power supply unit
Plug-in power supply / adapter for countries - optional	EURO

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	350×240×160 mm
Shipping method	Parcel service
Net weight approx.	1,6 kg
Gross weight approx.	1,8 kg
Shipping weight	2,7 kg
Shipping volume weight	2,7 kg

Product information

GTIN/EAN number	4045761478336
Trade number	90119000
REACH classification	no prohibited substances