

KERN OBB-A 1672

Filter geel

KERN



Pictogrammen

STANDARD



Categorie

Merk	Optics
Productcategorie	Optische Instrument
Productgroep	Microscoop filter

Ontwerp

Filter - type	Standard
Filter - kleur	geel
Filter - diameter	23 mm

Milieuomstandigheden

Opslagtemperatuur [Min]	-5 °C
Opslagtemperatuur [Max]	40 °C

Verpakking & verzending

Levertijden	1 d
Afmetingen verpakking (B×D×H)	70×45×45 mm
Verzendmethode	Pakketdienst
Nettogewicht ca.	0,05 kg
Brutogewicht ca.	0,05 kg
Verzendgewicht	0,1 kg
Verzendvolume (gewicht)	0,1 kg

Productinformatie

GTIN/EAN-nummer	4045761494459
Goederencode	85044095
REACH classificatie	geen verboden stoffen