

# KERN ORA-A2 107

Leather sleeve for analogue refractometers, black

**KERN**



## Category

|                  |                    |
|------------------|--------------------|
| Brand            | Optics             |
| Product category | Optical Instrument |
| Product group    | Leather sleeve     |

## Environmental conditions

|                           |       |
|---------------------------|-------|
| Storage temperature [Min] | -5 °C |
| Storage temperature [Max] | 40 °C |

## Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Readability force [d] (N)    | 1 d            |
| Dimensions packaging (W×D×H) | 200×150×70 mm  |
| Net weight                   | 0,15 kg        |
| Shipping method              | Parcel service |
| Net weight approx.           | 0,15 kg        |
| Gross weight approx.         | 0,30 kg        |
| Shipping weight              | 0,4 kg         |

## Pictograms

### STANDARD

