

# KERN OZB-A4204

Conversion objective lens 1,5 x

**KERN**



## Category

Brand	Optics
Product category	Optical Instrument
Product group	Conversion objective lens

## Construction

Anti fungus	✓
-------------	---

## Objective

Objective magnification	1,5 x
-------------------------	-------

## Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	55×55×25 mm
Net weight	0,7 kg
Shipping method	Parcel service
Net weight approx.	0,70 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg

## Pictograms

STANDARD

