

# KERN OZB-A45 15

# KERN

Double gooseneck illumination for best illumination of your sample



## Pictograms

### STANDARD



## Category

Brand	Optics
Product category	Microscope
Product group	Microscope illumination
Product family	OZB-IF

## Zulassung

CE mark	✓
---------	---

## Construction

Dimension (W×D×H)	120×370×250 mm
Beleuchtungsart	LED
Illuminance	6 W
Colour temperature	5600-6300 K
Gooseneck length	300 mm

## Power Supply

Supplied power supply	Power supply unit
Plug-in power supply / adapter for countries - optional	UK

## Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	300×250×250 mm
Net weight	1,74867 kg
Shipping method	Parcel service
Net weight approx.	1,8 kg
Gross weight approx.	2,0 kg
Shipping weight	3,8 kg