

# KERN OZB-A45 16

# KERN

Double gooseneck illumination for best illumination of your sample



Gross weight approx. 4,0 kg  
Shipping weight 4 kg

## Pictograms

## STANDARD



## Category

Brand	Optics
Product category	Microscope
Product group	Microscope illumination
Product family	OZB-IF

## Approval

CE mark	✓
---------	---

## Construction

Dimension (W×D×H)	100×200×150 mm
Beleuchtungsart	LED
Illuminance	20 W
Colour temperature	6400 K
Illumination dimmable	✓
Gooseneck length	540 mm

## Power Supply

Supplied power supply	Power supply unit
Plug-in power supply / adapter for countries - optional	EURO AUS UK US CH

## Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	300×250×250 mm
Net weight	3,64517 kg
Shipping method	Parcel service
Net weight approx.	3,8 kg