

# KERN OZB-A453 1

Gooseneck Illumination, 490 mm

**KERN**



## Category

Brand	Optics
Product category	Optical Instrument
Product group	Microscope illumination

## Construction

Gooseneck length	490 mm
Beleuchtungsart	Halogen

## Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	280×185×50 mm
Shipping method	Parcel service
Net weight approx.	0,70 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg

## Pictograms

### STANDARD

