


# KERN OZB-A4605

**KERN**

Mechanical specimen stage dimensions WxD 180×155 mm, path: 100x86 mm, (incident) + (transmitted light)



## Category

Brand	Optics
Product category	Optical Instrument
Product group	Microscope stage

## Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	260×240×100 mm
Net weight	1,5 kg
Shipping method	Parcel service
Net weight approx.	1,6 kg
Gross weight approx.	2,0 kg
Shipping weight	2 kg

## Pictograms

### STANDARD

