

KERN OZB-A4643

KERN

Conversion objective 2x

Pictograms

STANDARD



Category

Brand	Optics
Product category	Optical Instrument
Product group	Conversion objective

Construction

Anti fungus	✓
-------------	---

Objective

Objective magnification	2x
-------------------------	----

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	55×55×25 mm
Shipping method	Parcel service
Net weight approx.	0,70 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg

Product information

GTIN/EAN number	4045761191037
Trade number	90119000
REACH classification	no prohibited substances