

KERN OZB-A52 1 1

KERN

Individuality, variety and flexible working through our modular construction system ▶ PREMIUM-Universal stands



Pictograms

STANDARD



Category

Brand	Optics
Product category	Microscope
Product group	Stereo Microscope- Stand
Product family	OZB-UP

Construction

Dimension (W×D×H)	80×614×515 mm
-------------------	---------------

Stereomicroscope stand

Stand type	universal (telescopic arm with clamp)
Diameter mounting universal stand (Ø)	37 mm
Diameter mounting microscope head (Ø)	76 mm
Arm length	614 mm

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	680×330×260 mm
Net weight	5,8472 kg
Shipping method	Parcel service
Net weight approx.	6 kg
Gross weight approx.	7 kg
Shipping weight	11,7 kg