

# KERN OZB-A56 13

Achromatic conversion objective lens 0,7 x

# KERN



## Pictograms

### STANDARD



## Category

Brand	Optics
Product category	Optical Instrument
Product group	Conversion objective lens

## Construction

Anti fungus	✓
-------------	---

## Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	55×55×25 mm
Net weight	0,5 kg
Shipping method	Parcel service
Net weight approx.	0,50 kg
Gross weight approx.	0,70 kg
Shipping weight	0,7 kg

## Objective

Objective quality	Achromatic
Objective magnification	0,7 x
Objective working distance	145 mm