

KERN OZB-A56 17

KERN

Achromatic objective, 0,5x / W.D. (183 mm)



Trade number

90119000

REACH classification

no prohibited substances

Pictograms

STANDARD



Date of creation: 03.10.2026

Category

Brand	Optics
Product category	Optical Instrument
Product group	Microscope objective achromatic

Construction

Anti fungus	✓
-------------	---

Objective

Objective magnification	0,5x
Objective quality	Achromatic
Objective working distance	183 mm

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	55×55×25 mm
Shipping method	Parcel service
Net weight approx.	0,50 kg
Gross weight approx.	0,70 kg
Shipping weight	0,7 kg
Shipping volume weight	0,1 kg

Product information

GTIN/EAN number	4045761504530
-----------------	---------------