

# KERN OZB-A6303

**KERN**

Individuality, variety and flexible working through our modular construction system ▶ ECO-Universal stands



## Pictograms

### STANDARD



## Category

Brand	Optics
Product category	Microscope
Product group	Stereo Microscope- Stand
Product family	OZB-UE

## Construction

Dimension (W×D×H)	770×85×300 mm
-------------------	---------------

## Stereomicroscope stand

Stand type	universal (spring loaded arm with clamp/holder/ coarse focusing knob)
------------	---

## Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	580×420×370 mm
Net weight	7 kg
Shipping method	Parcel service
Net weight approx.	7 kg
Gross weight approx.	9 kg
Shipping weight	18 kg