

KERN OZB-A7101EU

KERN

Professional illumination units guarantee outstanding, uniform and strong illumination



Gross weight approx.	1,00 kg
Shipping weight	1 kg

Pictograms

STANDARD



Category

Brand	Optics
Product category	Microscope
Product group	Microscope illumination
Product family	OZB-IR

Approval

CE mark	✓
---------	---

Construction

Dimension (W×D×H)	99×90×33 mm
Beleuchtungsart	LED
Illuminance	4,5 W
Colour temperature	6500-7000 K
Illumination dimmable	✓
Ring illumination diameter	6,2 cm

Power Supply

Supplied power supply	Power supply unit
Plug-in power supply / adapter for countries - optional	EURO

Environmental conditions

Storage temperature [Min]	0 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	250×150×80 mm
Net weight	0,9375 kg
Shipping method	Parcel service
Net weight approx.	0,95 kg