

# KERN PCJ 600-2M-Prepack

The complete solution for prepack processes, with optional verification



Software EasyTouch

SET-17  
PrePack

## Category

Brand	KERN
Product category	Balance software bundle
Product group	Prepack balance
Product family	PCJ-Prepack

## Measuring System

Weighing system	Strain gauge
Weighing capacity [Max]	600 g
Readability [d]	10 mg
Reproducibility	10 mg
Linearity	± 30 mg
USP Minimum weight of sample (k = 2, U = 0.1%)	20 g
Minimum weight [Min]	500 mg
Resolution	60.000
Tare range	600 g
Possible calibration points	0 g; 100 g; 200 g; 400 g; 600 g
Automatic adjustment - Standard temperature deviation	2 °C
Automatic adjustment - Standard time interval	8 h
Stabilization time	3 s
Warm-up time	120 min
Eccentric loading at 1/3 [Max]	30 mg
Conformity assessment possible	✓

## Construction

Dimension housing (W×D×H)	163×245×80 mm
Dimensions completely mounted (W×D×H)	163×245×80 mm
Dimensions weighing surface (W×D)	131×131 mm
Dimensions weighing surface	131×131 mm
Level indicator	✓
Levelling feet adjustable	✓

## Functions

Alibi memory	internal memory
Alibi memory - number of memory spaces	250.000
Counting resolution (laboratory conditions)	30.000
Counting reference weight can be entered	✓
Smallest piece weight when piece counting - laboratory conditions	20 mg
Smallest piece weight when piece counting - normal conditions	200 mg

## Interface

EasyTouch (ET) - Compatibility	Vollständig
--------------------------------	-------------

## Power Supply

Battery operating time	20 h
Battery capacity	220 mAh
Battery voltage	3 V
Rechargeable battery operating time - backlight on	24 h
Rechargeable battery operating time - backlight off	48 h
Rechargeable battery charging time	8 h

## Environmental conditions

Ambient temperature [Min]	15 °C
Ambient temperature [Max]	35 °C
Humidity of environment [Max]	80 %

## Approval

CE mark	✓
Conformity assessment possible ex works	✓
Verification scale interval [e]	100 mg

## Services (optional)

Article number for conformity assessment	965-216
Article number for DAkkS calibration	963-127
Article number for certificate of conformity (verification)	969-517

## Packing & Shipping

Delivery time	1 d
---------------	-----

# KERN PCJ 600-2M-Prepack




The complete solution for prepack processes, with optional verification

Dimensions packaging (WxDxH)	235×365×175 mm
Shipping method	Parcel service
Net weight approx.	2,0 kg
Gross weight approx.	2,6 kg
Shipping weight	3 kg

## Pictograms

### STANDARD



### OPTION



### FACTORY

