

KERN 307-03

Individual weight, OIML Class E1, 5 g, Knob, Stainless steel polished



Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	8×8×15 mm
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Optional cases

Case - Plastic article number	317-030-400
Case - Aluminium article number	317-030-600
Case - Wood article number	317-030-100

Pictograms

OPTION



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E1
Product family	307

Measuring System

Conformity assessment possible	✓
Test weight - volume (cm ³)	0,625
Test weight - volume uncertainty (cm ³)	0,0078
Test weight - Density (kg/m ³)	8000
Test weight - Density uncertainty (kg/m ³)	100
Test weight - Density [Min] (OIML)	7620
Test weight - Density [Max] (OIML)	8420
Test weight - nominal value	5 g
OIML-compliant	✓
Weight OIML class	E1
Permissible tolerance according to OIML (+/-)	0,016 mg

Construction

Shape weight	Knob
Dimensions (Ø×H)	8×15 mm
Material	Stainless steel polished

Services (optional)

Article number for DAkkS calibration	963-233
Article number for guaranteed plus-deviation	921-233