

Product brochure

08.08.2026

KERN 307-05

Individual weight, OIML Class E1, 20 g, Knob, stainless steel polished

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
315-243	Tweezers 315-243 (stainless steel 105 mm)
315-250	Tweezers 315-250 (stainless steel 150 mm)
317-050-100	Wood weight case 317-050-100 (for nominal values 20 g for classes E1+E2+F1 for design Cylindrical/compact)
317-050-400	Plastic weight case 317-050-400 (for nominal values 20 g for classes E1, E2 for design Cylindrical/compact)
317-050-600	Aluminium weight case 317-050-600 (for nominal values 20 g for classes E1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
921-235	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-235
962-235R	DAkKS Recalibration 962-235R
963-235	DAkKS calibration certificate 963-235

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days