

# KERN 307-06

Individual weight, OIML Class E1, 50 g, Knob, Stainless steel polished



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML E1
Product family	307

## Measuring System

Conformity assessment possible	✓
Test weight - volume (cm <sup>3</sup> )	6,25
Test weight - volume uncertainty (cm <sup>3</sup> )	0,0781
Test weight - Density (kg/m <sup>3</sup> )	8000
Test weight - Density uncertainty (kg/m <sup>3</sup> )	100
Test weight - Density [Min] (OIML)	7920
Test weight - Density [Max] (OIML)	8080
Test weight - nominal value	50 g
OIML-compliant	✓
Weight OIML class	E1
Permissible tolerance according to OIML (+/-)	0,03 mg

## Construction

Shape weight	Knob
Dimensions (Ø×H)	18×28 mm
Material	Stainless steel polished

## Services (optional)

Article number for DAkkS calibration	963-236
Article number for guaranteed plus-deviation	921-236

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	18×18×28 mm
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

## Optional cases

Case - Plastic article number	317-060-400
Case - Aluminium article number	317-060-600
Case - Wood article number	317-060-100

## Pictograms

### OPTION

