

Product brochure

08.08.2026

# KERN 307-07

Individual weight, OIML Class E1, 100 g, Knob, stainless steel polished

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

☎ +49 7433 9933-0  
✉ [info@kern-sohn.com](mailto:info@kern-sohn.com)  
🌐 [www.kern-sohn.com](http://www.kern-sohn.com)

## PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

## PICTOGRAMS

Standard

Option



## ACCESSORIES

| Model       | Description                                                                                                         |
|-------------|---------------------------------------------------------------------------------------------------------------------|
| 315-250     | Tweezers 315-250 (stainless steel   150 mm)                                                                         |
| 317-070-100 | Wood weight case 317-070-100 (for nominal values 100 g   for classes E1+E2+F1   for design Cylindrical/compact)     |
| 317-070-400 | Plastic weight case 317-070-400 (for nominal values 100 g   for classes E1, E2   for design Cylindrical/compact)    |
| 317-070-600 | Aluminium weight case 317-070-600 (for nominal values 100 g   for classes E1 - M3   for design Cylindrical/compact) |

## SERVICES

| Model    | Description                                                         |
|----------|---------------------------------------------------------------------|
| 920-350  | Marking by laser engraving (individual marking) KERN 920-350        |
| 920-351  | Marking of test weights (last digits of serial number) KERN 920-351 |
| 921-237  | Guaranteed Plus-deviation ( $> +1/3$ OIML) KERN 921-237             |
| 962-237R | DAkKS Recalibration 962-237R                                        |
| 963-237  | DAkKS calibration certificate 963-237                               |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days