

Product brochure

08.08.2026

KERN 308-22

Set of milligram weights, E1, 1 mg - 500 mg, Platelet, , in plastic case


ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
315-247	Tweezers 315-247 (stainless steel, plastic 136 mm)
318-270	Brush 318-270

SERVICES

Model	Description
962-250	DAkkS calibration certificate 962-250
962-250R	DAkkS Recalibration 962-250R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days