

KERN 308-42

KERN

Set of milligram weights, E1, 1 mg - 500 mg, Wire, Stainless steel polished, in plastic case



Optional cases

Case - Plastic article number 328-222-420

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	308

Measuring System

Test weight - nominal value	1 mg - 500 mg
OIML-compliant	✓
Weight OIML class	E1

Construction

Dimension (W×D×H)	136×53×22 mm
Shape weight	Wire
Material	Stainless steel polished
Number of weights	12
Case of weight set	in plastic case

Services (optional)

Article number for DAkkS calibration	962-250
--------------------------------------	---------

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	136×53×22 mm
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg