

Product brochure

08.08.2026

# KERN 313-056

Weight set, E2, 1 mg - 500 g, Knob, stainless steel polished, in aluminium case

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

☎ +49 7433 9933-0  
✉ [info@kern-sohn.com](mailto:info@kern-sohn.com)  
🌐 [www.kern-sohn.com](http://www.kern-sohn.com)



## FEATURES

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

## PICTOGRAMS

Standard

Option



## ACCESSORIES

| Model       | Description                                                                                                                    |
|-------------|--------------------------------------------------------------------------------------------------------------------------------|
| 313-050-600 | Aluminium weight case 313-050-600 (for nominal values 1 mg - 500 g   for classes E1 - M1   for design Cylindrical/compact)     |
| 315-242     | Tweezers 315-242 (polyamide   100 mm)                                                                                          |
| 315-243     | Tweezers 315-243 (stainless steel   105 mm)                                                                                    |
| 315-245     | Tweezers 315-245 (stainless steel   250 mm)                                                                                    |
| 315-247     | Tweezers 315-247 (stainless steel, plastic   136 mm)                                                                           |
| 315-250     | Tweezers 315-250 (stainless steel   150 mm)                                                                                    |
| 317-009-600 | Aluminium weight case 317-009-600 (for nominal values 1 mg - 500 mg (einzel)   for classes E1 - M3   for design Fracture gram) |
| 317-010-100 | Wood weight case 317-010-100 (for nominal values 1 g   for classes E1+E2+F1   for design Cylindrical/compact)                  |
| 317-010-600 | Aluminium weight case 317-010-600 (for nominal values 1 g   for classes E1 - M3   for design Cylindrical/compact)              |
| 317-020-100 | Wood weight case 317-020-100 (for nominal values 2 g   for classes E1+E2+F1   for design Cylindrical/compact)                  |
| 317-020-400 | Plastic weight case 317-020-400 (for nominal values 1 g - 2 g   for classes E1, E2   for design Cylindrical/compact)           |

| Model       | Description                                                                                                         |
|-------------|---------------------------------------------------------------------------------------------------------------------|
| 317-020-600 | Aluminium weight case 317-020-600 (for nominal values 2 g   for classes E1 - M3   for design Cylindrical/compact)   |
| 317-030-100 | Wood weight case 317-030-100 (for nominal values 5 g   for classes E1+E2+F1   for design Cylindrical/compact)       |
| 317-030-400 | Plastic weight case 317-030-400 (for nominal values 5 g   for classes E1, E2   for design Cylindrical/compact)      |
| 317-030-600 | Aluminium weight case 317-030-600 (for nominal values 5 g   for classes E1 - M3   for design Cylindrical/compact)   |
| 317-040-100 | Wood weight case 317-040-100 (for nominal values 10 g   for classes E1+E2+F1   for design Cylindrical/compact)      |
| 317-040-400 | Plastic weight case 317-040-400 (for nominal values 10 g   for classes E1, E2   for design Cylindrical/compact)     |
| 317-040-600 | Aluminium weight case 317-040-600 (for nominal values 10 g   for classes E1 - M3   for design Cylindrical/compact)  |
| 317-050-100 | Wood weight case 317-050-100 (for nominal values 20 g   for classes E1+E2+F1   for design Cylindrical/compact)      |
| 317-050-400 | Plastic weight case 317-050-400 (for nominal values 20 g   for classes E1, E2   for design Cylindrical/compact)     |
| 317-050-600 | Aluminium weight case 317-050-600 (for nominal values 20 g   for classes E1 - M3   for design Cylindrical/compact)  |
| 317-060-100 | Wood weight case 317-060-100 (for nominal values 50 g   for classes E1+E2+F1   for design Cylindrical/compact)      |
| 317-060-400 | Plastic weight case 317-060-400 (for nominal values 50 g   for classes E1, E2   for design Cylindrical/compact)     |
| 317-060-600 | Aluminium weight case 317-060-600 (for nominal values 50 g   for classes E1 - M3   for design Cylindrical/compact)  |
| 317-070-100 | Wood weight case 317-070-100 (for nominal values 100 g   for classes E1+E2+F1   for design Cylindrical/compact)     |
| 317-070-400 | Plastic weight case 317-070-400 (for nominal values 100 g   for classes E1, E2   for design Cylindrical/compact)    |
| 317-070-600 | Aluminium weight case 317-070-600 (for nominal values 100 g   for classes E1 - M3   for design Cylindrical/compact) |
| 317-080-100 | Wood weight case 317-080-100 (for nominal values 200 g   for classes E1+E2+F1   for design Cylindrical/compact)     |
| 317-080-400 | Plastic weight case 317-080-400 (for nominal values 200 g   for classes E1, E2   for design Cylindrical/compact)    |
| 317-080-600 | Aluminium weight case 317-080-600 (for nominal values 200 g   for classes E1 - M3   for design Cylindrical/compact) |

| Model       | Description                                                                                                                        |
|-------------|------------------------------------------------------------------------------------------------------------------------------------|
| 317-090-100 | Wood weight case 317-090-100 (for nominal values 500 g   for classes E1+E2+F1   for design Cylindrical/compact)                    |
| 317-090-400 | Plastic weight case 317-090-400 (for nominal values 500 g   for classes E1, E2   for design Cylindrical/compact)                   |
| 317-090-600 | Aluminium weight case 317-090-600 (for nominal values 500 g   for classes E1 - M3   for design Cylindrical/compact)                |
| 317-280     | Gloves 317-280 (cotton, 1 pair)                                                                                                    |
| 318-01      | Test weight 318-01 (OIML Class E2, Nominal value 1 mg, Platelet, stainless steel polished)                                         |
| 318-02      | Test weight 318-02 (OIML Class E2, Nominal value 2 mg, Platelet, stainless steel polished)                                         |
| 318-03      | Test weight 318-03 (OIML Class E2, Nominal value 5 mg, Platelet, stainless steel polished)                                         |
| 318-04      | Test weight 318-04 (OIML Class E2, Nominal value 10 mg, Platelet, stainless steel polished)                                        |
| 318-05      | Test weight 318-05 (OIML Class E2, Nominal value 20 mg, Platelet, stainless steel polished)                                        |
| 318-06      | Test weight 318-06 (OIML Class E2, Nominal value 50 mg, Platelet, stainless steel polished)                                        |
| 318-07      | Test weight 318-07 (OIML Class E2, Nominal value 100 mg, Platelet, stainless steel polished)                                       |
| 318-08      | Test weight 318-08 (OIML Class E2, Nominal value 200 mg, Platelet, stainless steel polished)                                       |
| 318-09      | Test weight 318-09 (OIML Class E2, Nominal value 500 mg, Platelet, stainless steel polished)                                       |
| 318-270     | Brush 318-270                                                                                                                      |
| 328-222-420 | Weight case 328-222-420 (for nominal values 1 mg - 5 g   for classes E1 - M2   for design Cylindrical/compact)                     |
| 338-090-200 | Wood weight case 338-090-200 (for nominal values 1 mg - 500 mg (einzel)   for classes E1 - M2   for design Cylindrical/compact)    |
| 347-009-400 | Plastic weight case 347-009-400 (for nominal values 1 mg - 500 mg (einzel)   for classes E1 - M2   for design Cylindrical/compact) |

## SERVICES

| Model    | Description                                         |
|----------|-----------------------------------------------------|
| 920-705  | Marking (Last digits of serial number) KERN 920-705 |
| 962-304  | DAkkS calibration certificate 962-304               |
| 962-304R | DAkkS Recalibration 962-304R                        |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days